

# Coloured Caps for Cryovials

**Maxxline Coloured Caps for cryovials are supplied in a wide array of vibrant colours designed to make it effortless to identify the vials with various samples.**

These colourful inserts are the perfect complement to Maxxline's range of cryogenic vials, helping you to optimise your cryopreservation processes.

## Flexibility in Sample Processing


Maxxline's Coloured Caps come in a rainbow of colours, enabling customers to use them in a wide variety of sample sorting algorithms. Colour options include red, blue, yellow, green and white. These caps can be ordered as single-colour or multicolour packs with a mixture of up to five assorted colours.

## Durable Cap Inserts

Maxxline Coloured Caps are made from the same high-quality blend of polypropylene as the cryogenic vials and designed to fit snugly in the caps of the vials without falling out during use. They can withstand cryopreservation temperatures as low as -196°C.

## Cryo Vial Cap Inserts Free from Biological Contamination

Just like Maxxline Cryogenic Vials, the Coloured Caps are certified free from common biological contaminants such as endotoxins, DNA, RNases/DNases, and PCR inhibitors, among others. This protects the precious samples stored in the vials from unwanted contamination that might degrade sample quality.



## Ordering Information

Cat. No.	Description	Colour	Packaging
MLTR	Coloured Caps	Red	500 pcs.
MLTB	Coloured Caps	Blue	500 pcs.
MLTY	Coloured Caps	Yellow	500 pcs.
MLTG	Coloured Caps	Green	500 pcs.
MLTW	Coloured Caps	White	500 pcs.
MLTPAC	Coloured Caps	5 different colours	500 pcs.
MLD-50	Cryovial decapping tool		50 pcs.

