

Cardboard Cryogenic Boxes 136x136mm

Maxxline's 136x136mm Cardboard Cryogenic Boxes offer extended capacity for cryopreservation storage. The large 136x136mm interior extends the storage capacity of the 133x133mm Cryogenic Boxes, making them an excellent choice for laboratories with larger sample storage demands.

Temperature resistant Cryogenic Boxes

Our unique manufacturing process ensures that Maxxline Cryogenic Boxes hold up to the exacting demands of cryovial storage. Each box is meticulously crafted to withstand temperatures as low as -80°C without loss of integrity. Maxxline Cardboard Boxes have a water-repellent coating that reduces any moisture-induced sogginess.

Adaptable Boxes

As Maxxline Cardboard Boxes are supplied with or without dividers, they can be used to store a wide range of laboratory items. From routine biobanking to cryopreservation of ovarian tissue, these boxes are great for a variety of freezer applications.



Supplied in heights from 32mm to 130mm, they accommodate variable specimen tube heights.

Optional Cryo Vial Box dividers

Maxxline Cardboard Cryogenic Boxes can be partitioned with a range of optional grid dividers (7x7 to 16x16) making it easy to sort samples in each box and increasing flexibility for the user.

Highlights:

- Eco-friendly Cryogenic Storage Boxes with a low carbon footprint
- Cardboard Boxes with resistance to alcohols and mild organic solvents
- Thermostable and water-repellent
- Each box can be fitted with optional grid inserts for more efficient use of space
- Five box height variants fit a broad range of tubes (32mm, 50mm, 75mm, 100mm, and 130mm)
- Five colour variants for visual sample cataloguing (white, blue, green, red, and yellow)



Ordering Information

Bestseller Cardboard Cryogenic Boxes 136x136mm


Cat. No.	Description	Dimensions	Sales Unit
ML22500	Cardboard box, 50mm, white	136x136x50mm	36 pcs.
ML22684	Dividers for 50mm boxes, grid 9x9	136x136x30mm	36 pcs.
ML22683	Dividers for 50mm boxes, grid 10x10	136x136x30mm	36 pcs.

Cardboard Cryogenic Boxes 136x136mm with Dividers

Cat. No.	Description	Colour	Dimensions	Sales Unit
ML22231	Cardboard box B9 incl. 9x9 cells	White	136x136x50mm	36 pcs.
ML22233	Cardboard box B10 incl. 10x10 cells	White	136x136x50mm	36 pcs.
ML22933	Cardboard box w/print incl. 9x9 dividers	White	136x136x50mm	36 pcs.

Cardboard Cryogenic Boxes 136x136mm without Dividers

Cat. No.	Description	Colour	Dimensions	Sales Unit
ML22490	Cardboard cryogenic box, 32mm	White	136x136x32mm	36 pcs.
ML22492	Cardboard cryogenic box, 32mm	Blue	136x136x32mm	36 pcs.
ML22494	Cardboard cryogenic box, 32mm	Red	136x136x32mm	36 pcs.
ML22496	Cardboard cryogenic box, 32mm	Green	136x136x32mm	36 pcs.
ML22498	Cardboard cryogenic box, 32mm	Yellow	136x136x32mm	36 pcs.
ML22500	Cardboard cryogenic box, 50mm	White	136x136x50mm	36 pcs.
ML22501	Cardboard cryogenic box, 50mm	Blue	136x136x50mm	36 pcs.
ML22502	Cardboard cryogenic box, 50mm	Red	136x136x50mm	36 pcs.
ML22503	Cardboard cryogenic box, 50mm	Green	136x136x50mm	36 pcs.
ML22504	Cardboard cryogenic box, 50mm	Yellow	136x136x50mm	36 pcs.
ML22525	Cardboard cryogenic box, 75mm	White	136x136x75mm	36 pcs.
ML22526	Cardboard cryogenic box, 75mm	Blue	136x136x75mm	36 pcs.
ML22527	Cardboard cryogenic box, 75mm	Red	136x136x75mm	36 pcs.
ML22528	Cardboard cryogenic box, 75mm	Green	136x136x75mm	36 pcs.
ML22529	Cardboard cryogenic box, 75mm	Yellow	136x136x75mm	36 pcs.
ML22540	Cardboard cryogenic box, 100mm	White	136x136x100mm	36 pcs.
ML22541	Cardboard cryogenic box, 100mm	Blue	136x136x100mm	36 pcs.
ML22542	Cardboard cryogenic box, 100mm	Red	136x136x100mm	36 pcs.
ML22543	Cardboard cryogenic box, 100mm	Green	136x136x100mm	36 pcs.
ML22544	Cardboard cryogenic box, 100mm	Yellow	136x136x100mm	36 pcs.
ML22560	Cardboard cryogenic box, 130mm	White	136x136x130mm	36 pcs.
ML22561	Cardboard cryogenic box, 130mm	Blue	136x136x130mm	36 pcs.
ML22562	Cardboard cryogenic box, 130mm	Red	136x136x130mm	36 pcs.
ML22563	Cardboard cryogenic box, 130mm	Green	136x136x130mm	36 pcs.
ML22564	Cardboard cryogenic box, 130mm	Yellow	136x136x130mm	36 pcs.



Dividers for Cardboard Cryogenic Boxes 133x133mm

Cat. No.	Description	Height	Sales Unit
ML22090	Dividers for 32mm boxes, grid 7x7 (vials < 17mm O.D.)	22mm	36 pcs.
ML22091	Dividers for 32mm boxes, grid 8x8 (vials < 15mm O.D.)	22mm	36 pcs.
ML22092	Dividers for 32mm boxes, grid 9x9 (vials < 13mm O.D.)	22mm	36 pcs.
ML22093	Dividers for 32mm boxes, grid 10x10 (vials < 12mm O.D.)	22mm	36 pcs.
ML22094	Dividers for 32mm boxes, grid 12x12 (vials < 10mm O.D.)	22mm	36 pcs.
ML22095	Dividers for 32mm boxes, grid 13x13 (vials < 9mm O.D.)	22mm	36 pcs.
ML22096	Dividers for 32mm boxes, grid 14x14 (vials < 8mm O.D.)	22mm	36 pcs.
ML22097	Dividers for 32mm boxes, grid 16x16 (vials < 7mm O.D.)	22mm	36 pcs.
ML22111	Dividers for 50mm boxes, grid 3x3 (vials < 41mm O.D.)	30mm	36 pcs.
ML22112	Dividers for 50mm boxes, grid 4x4 (vials < 31mm O.D.)	30mm	36 pcs.
ML22113	Dividers for 50mm boxes, grid 5x5 (vials < 24mm O.D.)	30mm	36 pcs.
ML22114	Dividers for 50mm boxes, grid 6x6 (vials < 20mm O.D.)	30mm	36 pcs.
ML22106	Dividers for 50mm boxes, grid 7x7 (vials < 17mm O.D.)	30mm	36 pcs.
ML22108	Dividers for 50mm boxes, grid 8x8 (vials < 15mm O.D.)	30mm	36 pcs.
ML22110	Dividers for 50mm boxes, grid 9x9 (vials < 13mm O.D.)	30mm	36 pcs.
ML22120	Dividers for 50mm boxes, grid 10x10 (vials < 12mm O.D.)	30mm	36 pcs.
ML22115	Dividers for 50mm boxes, grid 12x12 (vials < 10mm O.D.)	30mm	36 pcs.
ML22116	Dividers for 50mm boxes, grid 13x13 (vials < 9mm O.D.)	30mm	36 pcs.
ML22117	Dividers for 50mm boxes, grid 14x14 (vials < 8mm O.D.)	30mm	36 pcs.
ML22118	Dividers for 50mm boxes, grid 16x16 (vials < 7mm O.D.)	30mm	36 pcs.
ML22124	Dividers for 75mm boxes, grid 3x3 (vials < 41mm O.D.)	40mm	36 pcs.
ML22125	Dividers for 75mm boxes, grid 4x4 (vials < 31mm O.D.)	40mm	36 pcs.
ML22126	Dividers for 75mm boxes, grid 5x5 (vials < 24mm O.D.)	40mm	36 pcs.
ML22127	Dividers for 75mm boxes, grid 6x6 (vials < 20mm O.D.)	40mm	36 pcs.
ML22130	Dividers for 75mm boxes, grid 7x7 (vials < 17mm O.D.)	40mm	36 pcs.
ML22132	Dividers for 75mm boxes, grid 8x8 (vials < 15mm O.D.)	40mm	36 pcs.
ML22134	Dividers for 75mm boxes, grid 9x9 (vials < 13mm O.D.)	40mm	36 pcs.
ML22136	Dividers for 75mm boxes, grid 10x10 (vials < 12mm O.D.)	40mm	36 pcs.
ML22137	Dividers for 75mm boxes, grid 12x12 (vials < 10mm O.D.)	40mm	36 pcs.
ML22147	Dividers for 100 & 130mm boxes, grid 3x3 (vials < 41mm O.D.)	65mm	36 pcs.
ML22148	Dividers for 100 & 130mm boxes, grid 4x4 (vials < 31mm O.D.)	65mm	36 pcs.
ML22149	Dividers for 100 & 130mm boxes, grid 5x5 (vials < 24mm O.D.)	65mm	36 pcs.
ML22150	Dividers for 100 & 130mm boxes, grid 6x6 (vials < 20mm O.D.)	65mm	36 pcs.
ML22140	Dividers for 100 & 130mm boxes, grid 7x7 (vials < 17mm O.D.)	65mm	36 pcs.
ML22142	Dividers for 100 & 130mm boxes, grid 8x8 (vials < 15mm O.D.)	65mm	36 pcs.
ML22144	Dividers for 100 & 130mm boxes, grid 9x9 (vials < 13mm O.D.)	65mm	36 pcs.
ML22146	Dividers for 100 & 130mm boxes, grid 10x10 (vials < 12mm O.D.)	65mm	36 pcs.

