

Non-sterile **Cryovials** with Separate Lids

Maxxline non-sterile cryovials can be ordered with separate lids. Having the lids and cryovials in separate packaging would be an excellent option for high-throughput manual sample processing where sterile tubes are not required.

Just like other Maxxline cryovials, these non-sterile cryovials are manufactured in a class 8 and 9 cleanroom facility and are certified free of human DNA, DNases/RNases, and Pyrogens.

Ordering Information

Cat. No.	Description	Size	Sales Unit
MLC1NSP	Cryovial, non-sterile, bulk packed	1.0mL	1000 pcs. w/separate lid.
MLC2NSP	Cryovial, non-sterile, bulk packed	2.0mL	1000 pcs. w/separate lid
MLC4NSP	Cryovial, non-sterile, bulk packed	4.0mL	1000 pcs. w/separate lid
MLC5NSP	Cryovial, non-sterile, bulk packed	5.0mL	1000 pcs. w/separate lid

Non-sterile **Cryovials** with Separate Lids (without print)

Maxxline non-sterile cryovials without print are the most cost effective cryo tube product in our portfolio.

These cryovials have neither graduations nor printed writing area making them the perfect option for laboratories with their own labelling routines. The separately packaged cryovials and lids fit well into a manual sample handling process flow. With Maxxline non-sterile cryovials without print you can stretch your budget without compromising on quality.

Ordering Information

Cat. No.	Description	Size	Sales Unit
MLC05NS	Cryovial, non sterile - bulked packed w/o print	0.5mL	1000 pcs. w/separate lid
MLC1NS	Cryovial, non sterile - bulked packed w/o print	1.0mL	1000 pcs. w/separate lid
MLC2NS	Cryovial, non sterile - bulked packed w/o print	2.0mL	1000 pcs. w/separate lid
MLC4NS	Cryovial, non sterile - bulked packed w/o print	4.0mL	1000 pcs. w/separate lid
MLC5NS	Cryovial, non sterile - bulked packed w/o print	5.0mL	1000 pcs. w/separate lid