

# White Label PCR Consumables

The White Label PCR Consumables are available as regular profile tubes of 0.2mL volume capacity which is a great choice for standard PCR/qPCR thermocyclers, strips of 8-tubes with caps and 96-well plates.

All White Label PCR Consumables are made out of polypropylene, with a specially polished interior which significantly prevents any protein adhesion on the tube walls, enabling you to fully recover your sample.

Certified DNA, DNase, RNA, RNase and Pyrogen-free, White Label PCR Consumables are the perfect choice for the wide majority of applications.

Thanks to the high-quality polypropylene (PP) material from which they are constructed, the White Label PCR Consumables can also be autoclaved, if required.

## White Label PCR Consumables for the Ultimate Heat Distribution

Due to their thin-walled construction, the White Label PCR Consumables offer the ultimate transmission and distribution of heat. Designed to be universally compatible with the majority of the thermocyclers on the market, White Label PCR Consumables are the optimal choice for worldwide laboratory set-ups.

The V-shape bottom of the White Label PCR Consumables enables full sample recovery and along with the thin walls and leak-proof closure, the tube design minimizes any condensation.

The evaporation rate of the White Label PCR Tubes is less than 3%. In addition to these features, the ergonomic flat cap design enables one-hand opening and closing, increasing your work efficiency, comfort and throughput.

## Highlights:

- Thin-walled
- Even temperature transmission
- Minimal condensation
- Evaporation less than 3%
- Universal thermocycler compatibility
- Low and regular profile tubes
- One-hand operation
- Certified DNA, DNase, RNA, RNase and Pyrogen-free
- Autoclavable

## Ordering Information

Cat. No.	Description	Dimensions	Packaging
WL3010800	PCR tube 0,2mL, clear, flat cap	0,2mL	10x1000 pcs.
WL3037020	PCR strips of 8 tubes 0,2mL, attached flat cap	8x0,2mL	5x10x12 pcs.
WL3047010	(q)PCR plate 96-well, 0,2mL, non-skirted, clear	96x0,2mL	5x5x10 pcs.
WL3048020	(q)PCR plate sealing membrane (adhesive film)		10x10 pcs.

