

White Label Reagent Reservoirs

The White Label Reagent reservoirs are made of a unique blend of Polystyrene and provide a sturdy, universal reservoir solution to any global laboratory.

Large Sturdy Base with Graduation Marks

Available in three different sizes - 10mL, 25mL, 50mL - the White Label Reagent Reservoirs provide a large, sturdy base with the broad edges, standing firmly on your workbench.

The graduation marks provide you extended work comfort and enable accurate dispensing and sample volume identification.

Angled Interior and Spill-proof

The angled interior of the White Label Reagent Reservoirs supports the collection of the liquid and decreases reagent waste.

While processing liquid you don't need to worry about spillage because the reservoirs are spill-proof.

Compact Storage and Easy Pouring

To store the remaining liquid you can empty the reservoirs using one of the corners.

For compact storage you can stack the reservoirs easily and separate them by using the reservoir lip.

The White Label Reagent Reservoirs are ideal for using with all Single and Multichannel Pipettes.

To suit the various pipettes, the White Label Reagent Reservoirs come in three different sizes:

- 10mL, V-shaped bottom (ideal with 8- and 16-channel pipettes)
- 25mL, flat bottom (ideal with multichannel pipettes)
- 50mL, V-shaped bottom (ideal with 12-channel pipettes)

Highlights:

- Multiple volume options
- Large, sturdy base
- Spill-proof
- Perfect for multichannel pipettes
- Compact storage



Ordering Information

Cat. No.	Description	Dimensions	Packaging
WL3701119	Reagent Reservoirs 10mL, bulk	10mL	5x100 pcs.
WL3702119	Reagent Reservoirs 25mL, bulk	25mL	5x100 pcs.
WL3703119	Reagent Reservoirs 50mL, bulk	50mL	5x100 pcs.